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
**Catalog No.** CG128001, CG128005, CG128025

**1. Intended Use**

This SARS-CoV-2 & MP & RSV & ADV & Flu A/B Antigen Kit (Colloidal Gold) is a lateral flow rapid chromatographic immunoassay intended for the in vitro qualitative detection and differentiation of antigens from SARS-CoV-2, Mycoplasma Pneumoniae (MP), Respiratory Syncytial Virus (RSV), Adenovirus (ADV) and influenza A (Flu A) and influenza B (Flu B) in nasal, nasopharyngeal, and oropharyngeal (throat) swab specimens from individuals suspected of respiratory infections.

Results are for the identification of antigens of SARS-CoV-2, MP, RSV, ADV, influenza A and influenza B, but does not differentiate between SARS-CoV and SARS-CoV-2 viruses and is not intended to detect influenza C antigens. These antigens are generally detectable in upper respiratory specimens during the acute phase of infection. Positive results indicate the presence of viral or bacterial antigens, but correlation with patient history and other diagnostic information is necessary to determine the infection status. Positive results do not rule out infection or co-infection with other bacteria or viruses. The agent detected may not be the definite cause of disease.

Negative results should be treated as presumptive, do not rule out SARS-CoV-2, MP, RSV, ADV or influenza A/B infection, and should not be used as the sole basis for treatment or patient management decisions, including infection control decisions. It's recommended to confirm negative results with other assays, if necessary, for patient management.

For professional use only.

**2. Test Principle**

Colloidal gold based immunochromatography is applied. The test is a combination of SARS-CoV-2 antigen test, MP antigen test, RSV antigen test, ADV antigen test, and influenza A/B antigen test. Briefly, the SARS-CoV-2, MP, RSV, ADV or Flu A/B antigens in the sample bind with the colloidal gold (detector) labeled anti-SARS-CoV-2 / anti-MP / anti-RSV / anti-ADV / Flu A/B antibodies on the conjugate pad to form the immunocomplexes. When the complexes migrate to the test line (coated with anti-SARS-CoV-2 or anti-MP or anti-RSV, anti-ADV, anti-Flu A/B antibodies), the complexes that contain the SARS-CoV-2 or MP or RSV or ADV or Flu A/B antigens will be captured. Similarly, the colloidal gold labeled mouse antibody will be captured in the control lines (coated with goat anti-mouse IgG). The complexes containing SARS-CoV-2 or MP or RSV or ADV or Flu A/B antigen will produce a colored line in the specific test line region, indicating the sample is positive for the respective antigen. Similarly, a colored line will also appear in each of the control line region indicating the adequate volume of specimen has been added and the membrane wicking has occurred.

**3. Materials Provided**

Component	1T CG128001	5T CG128005	25T CG128025
Test Cartridge	1	5	25
Sampling Swab	1	5	25
Extraction Tube	1	5	25
Tube Holder	/	/	1
Instructions For Use	1	1	1

**4. Material Required but Not Supplied**

- 4.1 Timer
- 4.2 Any necessary personal protective equipment

**5. Storage and Stability**

Store the kit at 2 – 30°C in a dry place and avoid direct sunlight. Do not freeze. The unopened cartridges are stable until the expiry date printed on the labels. It's recommended to open the test cartridge pouch just before use. Once the pouch is opened, the test should be started within 60 minutes.

**6. Sample Collection**

6.1 Nasal, nasopharyngeal, and oropharyngeal (throat) swab samples are acceptable for testing with this kit.

6.2 **Proper sample collection and handling is critical to the performance of the kit.** When collecting the swab sample, use the swab supplied in the kit. Do NOT touch the swab tip when handling the swab sample.

(a). Nasal swab collection: while gently rotating, insert the swab into one nostril. The swab tip should be inserted up to 1 to 2 cm (or 1 cm for children aged 2-14 years) from the edge of the nostril. Rotate the swab 10 times against the nasal wall. Remove and repeat the sampling process using the same swab for the other nostril.

(b). Nasopharyngeal swab collection: insert swab through the nares perpendicular to the nose (face) until resistance is encountered and the fingers touch the nose. Leave the swab in place for 15 – 30 seconds. Rotate the swab 3 times and remove it from the nasopharynx.

(c). Oropharyngeal (throat) swab collection: insert swab into the posterior pharynx and tonsillar areas. Rub swab over both tonsillar pillars and posterior oropharynx and avoid touching the tongue, teeth, and gums.

6.3 After collection, the sample is recommended to be tested immediately. If this is not possible, the sample should be stored in a dry and sterilized container. The swab sample can be stored at 2-8°C for up to 8 hours or at -20°C for up to 7 days.

**7. Test Procedure**

7.1 Allow the kit to come to room temperature (15 – 30°C) before testing.

7.2 Remove the Test Cartridge from its package. Place the cartridge on a clean, flat and dry surface. Label the cartridge with patient ID.

7.3 Remove the foil from the top of the Extraction Tube.

7.4 Place and soak the patient swab into the Extraction Tube. Rotate the swab 10 times while squeezing the sides of the tube.

7.5 Leave the swab in the tube for 1 minute.

7.6 Stir to mix the contents well. Remove the swab while squeezing the sides of the tube.


7.7 Attach and tighten the dropper tip to the top of the tube. Swirl or tap the tube bottom to mix the contents well.

7.8 Invert the Extraction Tube, and then vertically add **3 drops** (around 75 µL) into each sample well of the cartridge.

7.9 Start the timer: leave the cartridge at room temperature for **15 minutes**.

7.10 After the 15-minute incubation, read the results. Do not read the results after 30 minutes (from addition of the sample).

Note: Make sure there is sufficient light when reading and interpreting the results. Results read before 15 minutes or after 30 minutes may be false negative, false positive, or invalid.

**8. Interpretation of Results**

For each of the 3 test lane, the Control (C) line must appear for the test to be valid regardless of the appearance of the test line(s). Failure to observe the control line indicates the results are not reliable. When this

occurs, check the operation procedure carefully, and test again with a new sample. If the problem recurs, contact your local distributor.

**(1) Positive Result**

If the Control (C) line and the test line(s) for any of COV or MP or RSV or ADV or FluA or FluB appear, the test is positive for the corresponding pathogen, i.e., COV, SARS-CoV-2; MP, Mycoplasma Pneumoniae, RSV, Respiratory Syncytial Virus; ADV, Adenovirus; FluA, Influenza A; FluB, Influenza B.

**Note:**

- The test lines (A, B or T line region) may vary in intensity depending on the amount of antigen detected. Any visible (even light or faint) T line(s) should be interpreted as positive.
- Co-infection with more than one pathogen is rare. If the results are positive for more than one antigen, e.g., influenza A, B or SARS-CoV-2, it's recommended to retest with a new sample.

**(2) Negative Result**

If the Control (C) lines are visible, but none of the test lines (COV, MP, RSV, ADV, FluA, FluB) appear, the test is negative.

**Note:**

- Negative results should be treated as presumptive. Negative results do not rule out infection with SARS-CoV-2, MP, RSV, ADV, Flu A or Flu B, and should not be used as the sole basis for treatment or patient management decisions.

**(3) Invalid Result**

If the Control (C) line is not visible, the test is invalid. Insufficient sample volume or incorrect operation are the most likely reasons for invalid results. Review the procedure and repeat the test with a new sample.

**9. Limitations**

- 9.1 The test is only for qualitative detection of SARS-CoV-2, MP, RSV, ADV, Flu A and Flu B antigens.
- 9.2 This test detects both viable and non-viable SARS-CoV-2, MP, RSV, ADV, Flu A and Flu B. Test performance depends on the amount of antigen in the sample and may or may not correlate with pathogen culture or molecular assay performed on the same sample.
- 9.3 A negative test result may occur if the level of antigen in a sample is below the detection limit of the test, if the sample was collected or handled improperly or if the pathogen has amino acid changes in the target epitope region recognized by the antibodies utilized in the test.
- 9.4 Failure to follow the Test Procedure may adversely affect test performance and/or invalidate the test result.
- 9.5 Test results obtained with the test must be used in conjunction with other clinical history, findings, and epidemiological data, etc.
- 9.6 Positive test results do not rule out co-infections with other pathogens, differentiate between SARS-CoV and SARS-CoV-2 or identify specific influenza subtypes.
- 9.7 Negative results do not rule out infection with SARS-CoV-2, MP, RSV, ADV, Flu A or Flu B, and should not be used as the sole basis for treatment or patient management decisions, including infection control decisions. Negative results must be combined with clinical observations, patient history and epidemiological information, and confirmed with a molecular assay, if necessary, for patient management.

**10. Warnings and Precautions**

- 10.1 For in vitro diagnostic use only.
- 10.2 If the package has been damaged, the label cannot be seen clearly or if the kit has expired, do not use the kit.
- 10.3 Read the instructions carefully before starting the test. To obtain accurate results, the instructions must be followed.
- 10.4 Extraction Buffer of different lots are not interchangeable. The results may not be reliable if reagents from different lots are mixed or used together.
- 10.5 The test cartridge is for single test and cannot be reused. Do not use expired cartridges.
- 10.6 Do not eat the desiccant.

- 10.7 Avoid exposure of the skin, eyes, mouth, mucous membranes to the Extraction Buffer. If the contact occurs, flush with large amounts of water.
- 10.8 Use of gloves during sampling and testing is recommended.
- 10.9 The samples, used reagents and consumables are potentially infectious waste and should be disposed of in accordance with national and local regulations.

**11. Symbols**

Consult instructions for use	Lot number	Do not reuse
Use-by date	In vitro diagnostic medical device	Temperature limit 2 - 30°C
Manufacturer	Tests per kit	Avoid sunshine
Date of manufacture	Catalog number	Authorized representative in the European Community
This product complies with the requirements of the European Directive 98/79/EC		

**Goldsite Diagnostics Inc.**

No. 1001, Building No.1, GOLDSITE Mansion,  
No. 9 Rongtian Road, Kengzi, Pingshan District,  
Shenzhen, 518122 China  
Tel: 86 755 26890807  
Fax: 86 755 26890799

**CMC MEDICAL DEVICES & DRUGS, S.L.**

C/ Horacio Lengo No 18, CP 29006, Málaga-Spain

IFU4542002, V1, 05/2024